

The Oxford® Partial Knee
The most widely used and clinically proven
partial knee system in the world

References

1. Kim, K.T. *et al.* A Prospective Analysis of Oxford Phase 3 Unicompartmental Knee Arthroplasty. *Orthopedics*. 30(5 Suppl): 15-18, 2007.
2. Svard, U. and Price, A. 20-year survival & 10 year clinical results of the Oxford medial UKA. Presentation. 73rd Annual AAOS Meeting. Chicago, IL. 2006.
3. Incavo, S.J. *et al.* Tibial Plateau Coverage in Total Knee Arthroplasty. *Clinical Orthopaedics and Related Research*. 299: 81-85, 1994.
4. Goodfellow, J.W. and O'Connor, J.J. The Mechanics of the Knee and Prosthesis Design. *Journal of Bone and Joint Surgery (Br)*. 60(3): 358-69, 1978.
5. Bankston, A.B. *et al.* Comparison of Polyethylene Wear in Machined Versus Molded Polyethylene. *Clinical Orthopaedics and Related Research*. 317: 37-43, 1995.
6. Clarke, I.C. *et al.* Hip Simulator Wear Testing ArCom® vs. Extruded Bar Polyethylene. Presentation. 7th Annual Conference on Techniques & Science for Successful Joint Arthroplasty. 1995.
7. Keys, G.W. Reduced Invasive Approach for Oxford II Medial Unicompartmental Knee Replacement - A Preliminary Study. *The Knee*. 6(3): 193-96, 1999.



Oxford® Partial Knee

THE
Oxford
Partial Knee

Biomet UK Ltd
Waterton Industrial Estate
Bridgend
South Wales
CF31 3XA
+44 [0] 1656 655221
+44 [0] 1656 645454

Disclaimer

Biomet UK Ltd., as the manufacturer of this device, does not practice medicine and does not recommend any particular surgical technique for use on a specific patient. The surgeon who performs any implant procedure is responsible for determining and utilising the appropriate techniques for implanting prosthesis in each individual patient.



THE OXFORD PARTIAL KNEE

- Less postoperative pain^{1*}
- Less bone removal^{1*}
- More natural motion^{1*}
- **Clinically Proven**
at 10, 15 and 20 years²

Tibial Component

- Anatomical shape for excellent bone coverage³



Femoral Component

- Conforming, spherical design minimizes contact stress throughout entire range of motion
- Curved inner geometry for minimal bone removal
- Multiple sizes for various patient anatomies

Mobile Meniscal Bearing

- Free-floating bearing designed to remain fully congruent with femoral component throughout entire range of motion⁴
- Only true mobile meniscal bearing knee system approved for use in the U.S.
- Increased wear resistance with ArCom[®] Direct Compression Molded polyethylene^{5,6}



Specifically Engineered

- Patented instrumentation allows for easy and accurate balancing of flexion and extension gaps
- Minimal incision to avoid quadriceps disruption and allow for more rapid recovery than TKA⁷
- Size specific instrumentation allows for bone removal in precise 1mm increments and accurate, reproducible results



¹*Compared to total knee replacement