



**Oxford Titanium Niobium Nitride**

Partial Knee Replacement

# The Oxford Partial Knee Replacement is now available with a Titanium Niobium Nitride coating.

- Components are coated by Physical Vapor Deposition (PVD) with a Titanium Niobium Nitride Coating
- Uses mobile bearing philosophy which reproduces natural kinematics<sup>1</sup>
- Utilises the standard Oxford instrumentation
- 5 sizes of femoral component to optimize patient fit
- 6 sizes of Alpha tibial component to allow better bone coverage
- Uses cemented Oxford component design, demonstrating clinical survivorship at 10, 15 and 20 years<sup>2</sup>
- Less bone removal<sup>1</sup>
- More natural motion<sup>1</sup>

## Femoral Implant Order Codes

Order Codes	Description
159530TNBN	Extra Small
154600TNBN	Small
154601TNBN	Medium
154602TNBN	Large
154603TNBN	Extra large

## Tibial Implant Order Codes

Order Codes	Description
154718TNBN	Size A, Left Medial
154719TNBN	Size A, Right Medial
154720TNBN	Size B, Left Medial
154721TNBN	Size B, Right Medial
154722TNBN	Size C, Left Medial
154723TNBN	Size C, Right Medial
154724TNBN	Size D, Left Medial
154725TNBN	Size D, Right Medial
154726TNBN	Size E, Left Medial
154727TNBN	Size E, Right Medial
154775TNBN	Size F, Left Medial
154776TNBN	Size F, Right Medial

## References

1. Kim, K.T. *et al.* A Prospective Analysis of Oxford Phase 3 Unicompartmental Knee Arthroplasty. *Orthopedics*. 30(5 Suppl): 15-18, 2007.
2. Svard, U. and Price, A. 20-year survival & 10 year clinical results of the Oxford medial UKA. Presentation. 73rd Annual AAOS Meeting. Chicago, IL. 2006.  
\* Compared to Total Knee Replacement

## Disclaimer

This publication and all content, artwork, photographs, names, logos and marks contained in it are protected by copyright, trademarks and other intellectual property rights owned by or licensed to Biomet UK Healthcare Ltd or its affiliates. This brochure is for information purposes and must not be used, copied or reproduced in whole or in part for any purposes other than marketing by Biomet UK Healthcare Ltd or its authorised representatives. Any other purposes are prohibited. For information about partial knee replacement, including indications, contraindications, warnings, precautions and potential adverse effects, see the implant specific package inserted on Biomet's website under Essential Patient Risk Information. Oxford is a registered trademark of Biomet UK Ltd.



**Legal Manufacturer:**  
Biomet U.K., Ltd.  
Waterton Industrial Estates,  
Bridgend, South Wales  
CF31 3XA,  
U.K.